Caution: Federal law restricts the sale of this device to by, or under the order of a dentist.

Product Description

CLEARX aligners are made from thermoformable polymer sheets that are used to incrementally move teeth using gentle force from a sequence of aligners to be worn by the patient over a designated period. The aligners are held in place by pressure and can be removed by the patient at any time.

CLEARX aligners contain a unique shape memory feature in which each aligner has both a default temporary configuration and a permanent configuration resulting in a reduction in the number of aligners required for tooth alignment treatment. A boosting step for each aligner is necessary to activate this component.

Material Composition

CLEARX aligners are composed of a three-layered thermoformable orthodontic material. The top and bottom layer is of cycloaliphatic polyester with a polyether polyurethane layer in between.

Indications for Use/Intended Use

CLEARX aligners are indicated for the alignment of teeth during orthodontic treatment of tooth malocclusion in patients with permanent dentition (i.e., all second molars).

Contraindications for Use

Aligners are contraindicated for use in patients with active periodontal disease, active caries, or temporomandibular disorder.

Warnings

Patients with a history of plastic allergies should not use this product. Damaged or broken aligners should not be used as small pieces of plastic can be swallowed or aspirated. Patients should inform their dentist immediately if an aligner is lost or damaged.

Adverse Reactions

In rare circumstances, an individual may be allergic to plastics. In that case remove aligners and discontinue treatment. Teeth may tend to return to their previous position after treatment, so a retainer may be necessary.

MRI Safety Information

CLEARX aligners are MR safe and electrically non-conductive, non-metallic, and nonmagnetic.

Instructions

- Each aligner in the CLEARX pouch consists of two treatment stages that should be used for two consecutive weeks. The aligner is received with the first stage treatment, while the second treatment stage is stored as the shape memory, which will only be activated after boosting.
- The patient should wear the aligner **directly as received** for the first week (without boosting).
- After wearing the aligner for one week, the patient is advised to boost

the aligner as described in the "How To Boost" section.

- After wearing **the boosted aligner** for one week, the patient need to progress to the next aligners pouch set and repeat the process. The patient is advised to never discard used aligners, they should be kept again in their pouches after use, as they may serve as templates or retainers in case of broken attachments or aligners.
- If an allergic reaction occurs while wearing aligners, discontinue use and report it to K Line immediately.

How to Boost

Boosting refers to the process where the aligner is heated to a specific triggering temperature, causing its shape to shift to the next aligner shape, known as the boosted form.

Booster refers to the activation device where the aligners are placed in with hot water to be activated.

Step 1: Preparing the Aligners

Carefully place the CLEARX aligners into the booster (open side up).

Step 2: Add Some Heat

Boil a kettle of water. Once boiled, immediately pour hot water into the booster up to the maximal fill line. Ensure that the CLEARX aligners are fully submerged.

Step 3: Give It Time to Activate

Close the lid of the CLEARX booster to maintain the optimal temperature and wait 10 minutes.

Step 4: Open Booster

Open the lid of the booster carefully and gently remove the CLEARX aligners with the tweezers.

Step 5: Cool Down

Rest the CLEARX aligners for a few minutes at room temperature to cool down. Once cooled, the aligners are ready to use.

Warning: The aligner should **NEVER** be microwaved or placed in boiling water on a stove top

Cleaning the Aligners

Aligners must NOT be cleaned with abrasive toothpaste or soaked in mouthwash or any other cleaning agent as this may scratch or dull the plastic. Aligners should be cleaned with a soft toothbrush and cool water to prevent warping.

Storage of Aligners

Aligners should not be stored near a heat source or in direct sunlight, as this can cause them to distort.

Inserting the Aligners

- Aligners must be rinsed before insertion into the mouth. Teeth should be brushed and flossed after every meal prior to inserting the aligners.
- It should be ensured that the correct aligner is used at each scheduled change by checking the aligner number marked on both the aligner and its pouch.
- Minor discrepancies between the teeth and a new aligner are normal since the teeth have not moved into their new position. If significant pain or discomfort is experienced, please check back to K Line support.

Removing the Aligners

- The aligner should be grasped with fingers starting at the back (the molars) and gently dislodged. When the aligner is released from the molars, it can be pulled away at the canine region.
- If resistance is felt, it may be necessary for the patient to repeat the previous step in several locations. Excessive force should not be used as it may damage the aligner.

Tooth Nerves

Excess pressure by the aligner can cause teeth nerve damage.

Reporting of Adverse Events

A space has been provided in the patient IFU for the practitioner to provide their contact information in case of an adverse event with the patient's aligners. Any serious incident or potential serious incident involving CLEARX aligners must be reported to K Line Europe GmbH.

HOW TO SUBMIT A CASE

Registration

Registration can be completed at <https://clearxportal.com/user/register>

If registration is successful, your account will be activated shortly. After that, the first case can already be submitted.

Submit Patient Information

Please submit the patient information, clinical photos, and X-rays via the upload section of our portal.

If you are not using an intraoral scanner, please send silicone impressions and bite registration, well padded, to:

**K Line Europe GmbH Bonner Straße 203 J-K 40589
Düsseldorf, Germany**

Extraoral Pictures

Please upload all pictures in the upload section of the portal.

1. Frontal with a smile
2. Frontal without a smile
3. Profile view

Intraoral Pictures

4. Upper jaw occlusal view
5. Frontal in occlusion
6. Lower jaw occlusal view

7. Slight open bite
8. Lateral view on the right
9. Lateral view on the left

Treatment Plan

Once all necessary information has been received, the treatment plan 3D simulation should be created within three to five business days. The treatment plan 3D simulation and comments will be found on the portal account.

Confirming the Case

After discussing the case with the patient, please confirm it on our portal. Any questions and comments can be added via the comment box.

After confirming the case it takes about 10-12 working days plus shipping until the aligners arrive at the shipping address.

ATTACHMENTS AND BUTTONS

Attachments Placement

1. Prepare the required instruments and materials (acid etch, bonding agent, composite filling, light curing device).
2. Take out the template from the ClearX kit delivered and fill out the indicated reservoirs for attachments with composite respectively.
3. Make sure that no composite is outside the reservoir.
4. Apply acid etch to enamel surface precisely where the attachments will be placed.

5. Apply bonding agent precisely on corresponding surfaces.
6. Use light curing to cure the applied bonding agent.
7. Place firmly and with care, this is to ensure the template fits well on the teeth and gums, especially on the indicated areas where attachments are required.
8. Light cure the composite filling while holding carefully the template in place.
9. Remove the template and clean off any remaining flashes and excess composite.

Buttons Placement

1. Prepare the required instruments and materials (acid etch, bonding agent, composite filling, light curing device, mini-mold and handle).
2. Apply acid etch gel to the mid cervical third of the tooth to be extruded. This can be altered for other movements indicated in the treatment plan.
3. Apply bonding agent to the indicated area where the button will be attached.
4. Light cure.
5. Place the mini-mold on the handle, make sure the wings of the mini-mold are horizontal to the handle itself.
6. Inject composite filling the mini-mold. The composite could be a flowable composite or brackets bonding composite.
7. Place handle and mini-mold correctly with the wings of the mini-mold

parallel to the teeth incisal edge.

8. Light cure while holding handle firmly towards tooth structure.
9. Remove handle and finish the button with diamond bur.

Retention

Retention is a constituent component of CLEARX treatment. After completion of the treatment, it stabilizes the teeth permanently in their desired position.

Why is Retention Important?

Immediately after treatment, the structures of the periodontium need time to regenerate and are initially unstable in the new position. The late growth of the jaw, which occurs beyond adolescence, or the breakthrough of wisdom teeth can also have a negative impact on the results achieved. The retainer prevents the teeth from moving back to their original position in these phases.

Types of Retention

There are two types of retention that can be used. In addition to the removable retainer there is a permanent lingual retainer. The permanent lingual retainer is fixed, and the duration will differ from case to case.

A removable retainer is part of our treatment and is provided for each case. Upon request, the permanent lingual retainer can also be obtained from us.

1. Removable Retainer

The removable retainer is made from a similar material as our ClearX aligners. The lifespan of the retainer depends on proper care. It should be

cleaned with cleaning tablets. When not in use, the patient should store it in the aligner box.

Wearing Time:

1 st - 6 th month	6 th - 9 th month	9 th - 12 th month	12 th - 36 th month
18 - 21 hours/day	8 - 12 hours/night	3 nights/week	1 night/week

2. Permanent Lingual Retainer

The permanent lingual retainer is attached at the back of the anterior teeth and is particularly beneficial for maintaining stability over a long period of time. From the outside, the permanent retainer is not visible. It is made of a medical metal wire or aesthetic fiberglass. After a short adaptation phase the permanent retainer does not disturb the patient while eating or talking. A permanent lingual retainer requires no active participation by the patient. It is firmly connected with the teeth and can be neither forgotten nor lost.

TROUBLESHOOTING

Check the Pressure:

After the boosting process, the patient should notice a distinct pressure or improved grip from the aligner on the teeth. If the aligner feels the same as it did before the boosting process, it's likely that the boosting process hasn't been successful and needs to be repeated.

Right Temperature:

Ensure the boiled water is poured into the booster immediately after boiling.

Waiting too long can reduce the temperature and affect the process.

Right Position:

Ensure the CLEARX aligners are placed into the booster open side up and fully submerged in the hot water.

Seal the Lid:

Always close the booster lid securely during the boosting process to maintain as much heat as possible.

Right Timing:

Wait for at least 10 minutes before removing the aligners.

If the patient followed all the steps correctly and still doesn't feel the pressure, it's essential to repeat the boosting process. Remember, the pressure or grip is a sign that your treatment is progressing as planned!

SYMBOL EXPLANATION

Symbol	Title	Description	Reference
	Manufacturer	Indicates the medical device manufacturer.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Use-by	Indicates the date after which the medical device is not to be used.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
LOT	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
REF	Reference number	Indicates the manufacturer's reference number so that the medical device can be identified.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Temperature limit	Indicates the temperature range at which the device need to be maintained.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements

Symbol	Title	Description	Reference
	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Prescription only	Requires a prescription in the United States of America Caution: In the United States of America, federal law restricts this device to sale or use by, or on the order of a dentist.	USA Code of Federal Regulations 21CFR Part 801 § 801.109 (b)(1)
	Magnetic resonance	CLEARX aligner is MR Safe.	FDA Guidance document: Testing and Labeling Medical Devices for Safety in the Magnetic Resonance (MR) Environment dated May 20, 2021.
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	EN ISO 15223 - 1:2021
	Non Sterile	Indicates a medical device that has not been subjected to a sterilization process.	EN ISO 15223 - 1:2021
	Medical device	Indicates the item is a medical device.	EN ISO 15223 - 1: 2021

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Düsseldorf, Germany ٤٠٥٨٩ | J-K ٢٠٣ K Line Europe GmbH| Bonner Straße